


Isolate the microscope from vibrations generated from the external environment



*Note:*  
AVT-1 is designed to be easy to use and almost maintenance free

**Technical Specifications AVT-1**

Overall dimensions (W x D x H)	1200 x 800 x 800 mm (47.25" x 31.5" x 31.5")
Weight	70 kg (154 lbs)
Material	Material Powder painted mild steel, Stainless steel and glass
Float dimensions (W x D x H)	540 x 340 mm (21.26" x 13.39")
Range (HZ) vibration could be eliminate while using this AVT	5,5-50 Hz



Designed in Denmark



Made in the E.U.